

ACD/LABS [ADVANCED CHEMISTRY DEVELOPMENT, INC.]

Accurate and Reproducible Quantitation of xC/UV/MS Data from Any Instrument

Anne Marie Smith
Baljit Bains

How Much Analyte is in the Sample?

Quantitative analysis addresses scientific challenges across various industries such as food, environment, pharmaceutical, and more. It is crucial to ensure the quality of a product by answering “*How much?*” This question is answered using essential quantitation xC/UV/MS workflows. Using MS in tandem with gas chromatography (GC) or liquid chromatography (LC) detects and quantitates analytes of interest in highly complex mixtures and allows for an efficient, selective, and sensitive separation.

In present-day multi-vendor laboratories, there are challenges in learning and navigating multiple applications with different interfaces and difficulty in standardizing data. Using software compatible with most leading vendor data formats minimizes limitations in instrument suitability and selection—enabling consistency and standardization of quantitative analysis.

MS Workbook Suite offers a fully integrated quantitative workflow where all xC/UV/MS data can be processed and quantitated in a single interface—ensuring consistently reliable and accurate results.

An Example Quantitation Workflow in MS Workbook Suite




Figure 1: The intuitive calibration and quantitation workflow in MS Workbook Suite.

Quantitation studies were designed within the software post-acquisition for a mixture prepared in triplicate (containing 4 compounds and 7 standards). Using the quantitation workflow to determine the amount of analyte present, 24 unknowns were quantified against the set of standards. The workflow consists of:

SET UP

Define sample information:

- *Sample type*—calibration, standard, blank, unknown sample
- *Type of trace*—extracted ion chromatogram (XIC), flat chromatogram (specified wavelength pre-acquisition), or DAD (user-specified wavelength during calibration)
- *Compound(s) for quantitation*—material A, B, C, or D

PROCESS

To determine the amount of analyte present, quantify the unknown samples against the set of standards. Define the processing parameters and identify the peaks for quantitation via retention time. Apply the peak detection and integration parameters to create a calibration curve and calculate the linear regression, standard deviation, and r^2 . Determine the suitability and validity of the statistical regression model by using the analysis of residuals (ANOVA table).




Figure 2. Chromatograms in MS Workbook Suite showing the associated structure and peaks for all 4 components, the internal standard, and the unknown at various concentrations.




Figure 3. Calibration curve for Material D in MS Workbook Suite, showing the regression line, standard deviation, residuals plot, amount of unknown (A3), and summary table (including r, r², equation, etc.)

REVIEW, STORE & REPORT

Automatically save the full set of data (raw and processed) to a database. It is possible to store and visualize data with its chemical context including chromatograms, calibration curves, and tabulated results. View or modify the results (i.e., alteration of processing parameters, or the addition of unknowns to be quantified) and track any changes made in the processing interface for data integrity purposes. Finally, you can create and share customized reports.

Processing and quantitating all xC/UV/MS data in MS Workbook Suite ensures connectivity between data and numerical results. Reproducible and accurate quantitative analysis allows confident determination of "how much" is in the sample ensuring the quality of final products whether food, water, or drugs.